

UNIVERSITÀ DEGLI STUDI DI MODENA E REGGIO EMILIA




A Cone Beam Computed Tomography Annotation Tool for Automatic Detection of the Inferior Alveolar Nerve Canal

Cristian Mercadante¹, Marco Cipriano¹, Federico Bolelli¹, Federico Pollastri¹, Mattia Di Bartolomeo², Alexandre Anesi³, and Costantino Grana¹

¹Department of Engineering "Enzo Ferrari", University of Modena and Reggio Emilia ²Surgery, Dentistry, Maternity and Infant Department, Unit of Dentistry and Maxillo-Facial Surgery, University of Verona ³Department of Medical and Surgical Sciences for Children & Adults, Cranio-maxillo-facial Surgery, University of Modena and Reggio Emilia

The Inferior Alveolar Nerve




- Knowing the position of the Inferior Alveolar Nerve (IAN) is essential in surgical operations.
- Traditional annotations are made on panoramic views using proprietary software.

A Novel Annotation Tool

The proposed graphic tool allows experts to load volumes from DICOM files, interactively evaluate the best cutting strategies, label the IAN from 2D slices and export them as 3D annotations.

🕼 IAN Anjanotation Tool - [Y:/work/datasets/maxillo/DICOM/anonymous/patient/DICOMD... – 🛛 🛛 🗙



• Labeling a 3D volume is tedious and discourages the creation of large datasets for further purposes (e.g. deep learning applications).

This work proposes a novel tool for creating 3D annotations of the IAN with minimum effort and high precision.


Cross-Sectional Views

Using the tool, an expert can:

- (a) Create and edit the panoramic base spline on dental arch;
- (b) Generate straight cutting planes orthogonal to the

Cross-Sectional Planes


- Tilted views can be obtained by following the slope of the canal.
- Those views are generated by tilting the cutting planes (1).




(1)


spline;

(c) Use planes to cut the volume and retrieve 2D CSVs. Annotate each view and get hints for the next one.



- The **slope** comes from a previous annotation on the panoramic view (2).
- Those images guarantee stable circular annotations even in critical positions (3).






Cut using straight planes vs tilted planes in the proximity of a tooth

3D Reprojection and saving

2D views can be re-projected to the original volume as **3D labels**.




Use case and results

- Annotation tool evaluation is performed generating three datasets from four patients and training the classic UNet3D architecture for a segmentation task.
- The network is fed with blocks of consecutive slices.

Traditional Annotated Volume


Our Annotation

Set of CSVs or 3D volumes can be exported in **different formats** for further purposes.



	2 Labels			3 Labels		
	GTD	GSD	ASD	GTD	GSD	ASD
P1 P2 P3 P4	$\begin{array}{c} 0.67 \\ 0.44 \\ 0.59 \\ 0.61 \end{array}$	$\begin{array}{c} 0.68 \\ 0.42 \\ 0.42 \\ 0.54 \end{array}$	$\begin{array}{c} 0.63 \\ 0.45 \\ 0.41 \\ 0.46 \end{array}$	$\begin{array}{c} 0.54 \\ 0.67 \\ 0.63 \\ 0.65 \end{array}$	$\begin{array}{c} 0.55 \\ 0.34 \\ 0.26 \\ 0.35 \end{array}$	$\begin{array}{c} 0.51 \\ 0.29 \\ 0.26 \\ 0.26 \end{array}$

Results of the segmentation task: IoU (IAN for 2 labels) and mIoU (internal IAN and contour for 3 labels) of our network trained with slices generated by guided tilted planes (GTD), guided straight planes (GSD), and unsupervised straight planes (ASD).



Qualitative result of a prediction. Original image on left, annotation from the tool in the middle, network prediction on the right.

 The annotation tool produces valid datasets for neural networks learning and generalization process.